

Cell Chemical Biology, Volume 26

Supplemental Information

Refactoring the Cryptic Streptopenazine

Biosynthetic Gene Cluster Unites Phenazine,

Polyketide, and Nonribosomal Peptide Biochemistry

Katherine D. Bauman, Jie Li, Kazuya Murata, Simone M. Mantovani, Samira Dahesh, Victor Nizet, Hanna Luhavaya, and Bradley S. Moore




Figure S1. Confirmation of the integrity of the captured *spz* BGC, related to Figures 1 and 2 and Table 1. **A.** Restriction digest of pSMM before and after transformation in *E. coli* using KpnI and KpnI + SpeI enzymes. KpnI digest, expected sizes (kb) = 23.9, 11.1, 8.2, 4.9, 4.1, 2.7, 1.9, 1.5, 1.1, 1.0, 0.5. KpnI + SpeI digest, expected sizes (kb) = 22.9, 8.9, 8.2, 4.9, 4.1, 2.7, 2.1, 1.9, 1.5, 1.1, 1.0, 0.5. **B.** Restriction digest of pKDB01, expected sizes (kb): 18.0, 13.7, 8.3, 4.9, 2.0, 1.6, and pKDB02, expected sizes (kb): 17.0, 13.3, 8.3, 4.9, 2.0, 1.5, with EcoRI + Stul. * indicates digest is blocked due to overlapping *dcm* methylation. Std – sequenced plasmid digested with EcoRI + Stul, used as a control, expected sizes (kb) = 10.3, 3.2, 1.3. Plasmid backbone represented by blue rectangle.





Figure S2. Engineering the *spz* BGC, related to Figures 2 and 3. A. Agarose gels of restriction digestions and corresponding maps of regulatory gene mutant plasmids. Enzyme pair *XhoI* + *NdeI* used for digestion. pKDB01 Δ *spz3*, expected sizes (kb): 10.6, 10.5, 8.7, 8.6, 3.2, 3.1, 1.9, 1.0, 0.1. pKDB02 Δ *spz3*, expected sizes (kb): 10.6, 10.5, 8.7, 8.6, 3.2, 2.1, 1.9, 1.0, 0.1. pKDB01-*ermE***p-spz24*, expected sizes (kb): 11.7, 10.5, 8.7, 5.2, 3.4, 3.2, 3.1, 1.9, 0.9, 0.2, 0.1. **B.** Restriction digest with *XhoI* + *NdeI* of refactored plasmids, agarose gel and corresponding plasmid map showing restriction sites: pKDB03, expected sizes (kb): 11.7, 10.5, 8.7, 7.6, 3.2, 3.1, 1.9, 1.1, 0.1; pKDB04, expected sizes (kb): 11.7, 10.5, 8.6, 7.5, 3.2, 3.1, 1.9, 1.1, 0.7, 0.3, 0.2, 0.1; pKDB05, expected sizes (kb): 10.9, 10.5, 8.7, 7.5, 3.2, 3.1, 1.9, 1.1, 0.9, 0.9, 0.7, 0.3, 0.3, 0.1. Plasmid backbone represented by blue rectangle.

A

Key, nodes observed in:

- M1146-pKDB03
- M1146-pKDB03, 04, 05


B

Figure S3. Comparative production of streptophenazines by the refactored *spz* BGC, related to Figures 2, 3, and 4. A. Streptophenazine metabolite cluster of nodes from the molecular network including MS/MS data of extracts of M1146-pKDB03, -pKDB04, and -pKDB05. Nodes are labeled with the corresponding precursor ion mass. **B.** Extracted ion chromatograms (m/z 411.2, 425.2, 439.2, 510.2), corresponding to isolated compounds in refactored pathway: (i) M1146-pKDB01, (ii) M1146-pKDB03, (iii) M1146-pKDB04, and (iv) M1146-pKDB05.





Figure S4. Molecular network showing production of streptopenazines, related to Figure 4. LC-MS/MS data of ethyl acetate extracts from CNB-091, M1146-pCAP03, M1146-pKDB01 (non-refactored *spz* BGC) and M1146-pKDB03 (refactored *spz* BGC) were used for network generation. Nodes corresponding to streptopenazines highlighted in red square.

A**B**

```

Spz7_KR      ...-D V Y W A H L P D T --V T P V E E I L A A L D D L V R A G K I L H A G L S N F P A W R V ...
Ery1_KR      ...---L G G I G D D V P L S A V F H A A A T L D D G T -V D T L T G E R I E R A S R A K V L ...
SlnA1_KR     ...---L L D R I P E A H P L T G V F H A A G V L D D G M -V G A L S A E R L D A V L R P K T D ...
RifA_KR      ...L E A V L R A I P A E H P L T A V I H T A G V L D D G V -V T E L T P D R L A T V R R P K V D ...
AmphJ_KR     ...---L L A S V P A E H P L T A V V H T A G V L D D G I -F P S L T P D R L D S V M R P K V D ...
TlmH_KR      ...---V L A Q I R S R G P I G G V V H A A G L L D D S I -L A N M T P E Q L H R V L R S K V D ...

```

Figure S5. Confirmation of stereochemistry at C-1' position of compound 18, related to Figure 1. A. Circular dichroism (CD) spectrum of streptopenazine G (**18**). **B.** Alignment of type I PKS KR domains and Spz7. Amino acid residues (LDD motif) defining Spz7 as type B KR are highlighted. Ery1 = erythromycin, SlnA1 = salinomycin, RifA = rifamycin, AmphJ = amphotericin, TlmH = thiolactomycin.




Figure S6. Bioinformatic analysis and gene deletion of *spz15*, related to Figures 5 and 6. **A.** Protein family (Pfam) analysis of *Spz15*. Identified CL00378 AMP-binding domain which is characteristic for ANL superfamily of enzymes that includes adenylation domains. **B.** Restriction digest of pKDB03 Δ *spz15* with *Xho*I and *Nde*I restriction enzymes. Expected sizes (kb): 11.7, 10.5, 7.5, 5.6, 3.2, 3.1, 2.3, 1.9, 1.1, 0.1. **C.** LC-MS chromatograms: (i) Base Peak Chromatogram (BPC) of M1146-pKDB03 Δ *spz15* *zoomed out 100x, (ii) Extracted Ion Chromatogram (EIC) (m/z 510.2, corresponding to compounds **20** and **21**) of M1146-pKDB03 Δ *spz15*, (iii) M1146-pKDB03 EIC (m/z 510.2), (iv) CNB-091 EIC (m/z 510.2), and (v) M1146-pCAP03 EIC (m/z 510.2).




Figure S7. antiSMASH predicted gene neighborhoods of *phz*-associated discrete adenylation proteins, related to Figure 5. Red gene = adenylation enzyme homolog found through blastp, orange genes = phenazine biosynthesis homologs. Spz15 = from *Streptomyces* sp. CNB-091, WP_097877165.1 = from *Streptomyces* sp. ms184, WP_086670704.1 = from *Streptomyces albovinaceus*, WP_073774768.1 = from *Streptomyces* sp. TSRI0445, WP_093779888.1 = from *Streptomyces* sp. yr375, WP_078880735.1 = from *Kitasatospora purpeofusca*, WP_063798176.1 = from *Streptomyces* sp. 150FB, WP_030586976.1 = from *Streptomyces anulatus*, WP_078871187.1 = from *Streptomyces caatingaensis*, WP_075016162.1 = from *Streptomyces rubidus*, WP_058941148.1 = from *Streptomyces kansasensis*.

Table S1, related to Figures 2, 3 and Method Details. Promoter cassette sequences used for refactoring. Bold = restriction site, italic = added for scar, blue = FRT site, red = antibiotic resistance, green = promoter, purple = *actII-ORF4*.

Cassette	Sequence
<i>sp44-p21</i>	<p> CTCGAGGGTGAACCGATCTCCTCGTTAGGGTCAACCCAGACTTTACAACCCGCACAGCATGT TGTCAAAGCAGAGACGGTTTGAATGTGAACAGCCACTATCATATGTGCAGTTCGAAGTTCCTAT TCTCTAGAAAGTATAGGAACTTCGGTTCATGTGCAGCTCCATCAGCAAAGGGGATGATAAGTT TATCAACCACCGACTATTTGCAACAGTGCCGTTGATCGTGCTATGATCGACTGATGTCATCAGCG GTGGAGTGCAATGTCGTGCAATACGAATGGCGAAAAGCCGAGCTCATCGGTCAGCTTCTCAAC CTTGGGGTTACCCCCGGCGGTGTGCTGCTGGTCCACAGCTCCTTCCGTAGCGTCCGGCCCT CGAAGATGGGCCACTTGGACTGATCGAGGCCCTGCGTGCTGCGCTGGGTCCGGGAGGGACG CTCGTCATGCCCTCGTGGTCAGGTCTGGACGACGAGCCGTTGATCCTGCCACGTCCGCCGT TACACCGGACCTTGAGTTGTCTTGACACATTCTGGCGCCTGCCAAATGTAAAGCGCAGCGC CCATCCATTTGCCTTGGCGGACGCGGGCCACAGGCAGAGCAGATCATCTCTGATCCATTGCC CCTGCCACCTCACTCGCCTGCAAGCCCGGTGCCCGTGTCCATGACTCGATGGGCAGGTA TTCTCCTCGGCGTGGGACACGATGCCAACACGACGCTGCATCTTGCCGAGTTGATGGCAAAG GTTCCCTATGGGGTCCGAGACACTGCACCATTCTCAGGATGGCAAGTTGGTACGCGTCGAT TATCTCGAGAATGACCACTGCTGTGAGCGCTTGCCTTGGCGGACAGGTGGCTCAAGGAGAA GAGCCTTCAGAAGGAAGGTCCAGTCGGTCATGCCTTGTCTCGGTTGATCCGCTCCCGCGACAT TGTGGCGACAGCCCTGGGTCAACTGGGCCGAGATCCGTTGATCTTCTGCATCCGCCAGAGG CGGGATGCGAAGAATGCGATGCCGCTCGCCAGTCGATTGGCTGAGGAAGTTCCTATTCTCTAGA AAGTATAGGAACTTCAAGCTTTGCTCGAGTGTGCGGGCTCTAACACGTCCTAGTATGGTAGGA TGAGCAATCTAGTCGAGCAACGGAGGTACGGACCATATG </p>
<i>actIp</i>	<p> CTCGAGGGTTCATGTGCAGCTCCATCAGCAAAGGGGATGATAAGTTTATCAACCACCGACTAT TTGCAACAGTGCCGTTGATCGTGCTATGATCGACTGATGTCATCAGCGGTGGAGTGCAATGTC GTGCAATACGAATGGCGAAAAGCCGAGCTCATCGGTCAGCTTCTCAACCTTGGGGTTACCCCC GGCGGTGTGCTGCTGGTCCACAGCTCCTTCCGTAGCGTCCGGCCCTCGAAGATGGGCCACT TGGACTGATCGAGGCCCTGCGTGCTGCGCTGGGTCCGGGAGGGACGCTCGTCAATGCCCTCG TGGTCAGGTCTGGACGACGAGCCGTTGATCCTGCCACGTCCGCCGTTACACCGGACCTTGG AGTTGTCTTGACACATTCTGGCGCCTGCCAAATGTAAAGCGCAGCGCCCATCCATTTGCCTT GCAGCGGGGGCCACAGGCAGAGCAGATCATCTGATCCATTGCCCTGCCACCTCACTC GCCTGCAAGCCCGGTCCGCTGTCATGAACTCGATGGGCAGGTACTTCTCCTCGGCGTGG GACACGATGCCAACACGACGCTGCATCTTGCCGAGTTGATGGCAAAGGTTCCCTATGGGGTG CCGAGACACTGCACCATTCTCAGGATGGCAAGTTGGTACGCGTCGATTATCTCGAGAATGAC CACTGCTGTGAGCGCTTGCCTTGGCGGACAGGTGGCTCAAGGAGAAGAGCCTTCAGAAGGA AGGTCCAGTCGGTCATGCCTTGTCTCGGTTGATCCGCTCCCGCGACATTGTGGCGACAGCCCT GGGTCAACTGGGCCGAGATCCGTTGATCTTCTGCATCCGCCAGAGGCGGGATGCGAAGAAT GATGCCGCTCGCCAGTCGATTGGCTGACCATATGCCACTGCCTCTCGGTAATAATCCAGCAAA AATTAATCAGTGACGCTCGCTGCACTGATTAATTTTTGATCAATAGGAGATCGCTTGTGACGGC AAGCACATTGAAATCTGTTGAGTAGGCCTGTTATTGTCCGCCCCAGGAGACGGAGAATCTCGA CGGGGGCGCAGATGAGATTCAACTTATTGGGACGTGTCCATGTAATCACCGATGCGGGATGTG TAATTCGGCTTAAATCCTCGAAGGCGACCCAGCTCCTGGTGCTGCTGCTCCTCAGGCGGCACG AGGTGGTGGGATCGGGGGTGTCTATCGAGGAGTTGTGGGCGGACCACCCGCCCGCAGCGC CATGACGACGCTGCAGACGTACGTGTACCACACCCGCCGGCTGCTGGGGGAGCACCGGGTG ACGAGCGACGACCGGGAATTGGTCTGACCCAGCCGCGCCGGCTACTTCGCCCTGATCGACGA GGACGAACTCGACGTCGCGGTGCGCCGAGCGTGTATCCGCACCGGCGCCGGCTGCTCGAG GAGAACCGGCTCGAGGAGGCGCTCGCCTCGTTGGACGCGGGACTGGATCTCTGGCGAGGCC CGGCGCTGTCCACCGTACCGTGCGGCGGGTGTCTGAAAGCAATATCGCGCACCTGGAAGAG CTGCGGCTTTTTGGAATGCAGCTCCGTATCGACGCGAATTGGCGGCTGGGCAGAATAGGGCC GATGATTCGGAACTCCGGTCCCTGGTAATTTTCGCATCCGCTGAACGAGACCCTGCACGCCAA ACTGATGGGCGCGCTCTGTGAGATGGGCAGGCGCGCCGAGGCGCTGGAATCGTATCGGAATC TCCGGCGGATACTGTCCGACGAAGTGGGGTGGATCCGACGCCGGAATCCAGCGTATGAC ATGGAATTTCAACGGTGAGAAGGTGCTCGTAGCACCGGTCCGTTGACCGGTTGAGCCGTTGAGCC CTATGTCTCTTAAGTGTTCCCCTCCCTGCCTCGTGGTCCCTCACGCGCTCAGCTTTGGGCGCC CGGCTCGAGCGGCGGTGCAAGGGAGATGGGGTCCGCTGGACGCGGGCGCCGGTGGATCCG GCATCGAGGGGTCCCGTATCGGCCTTCGAGCCTCCTTCGAGCCACGGGGCCGACGATGACGA CGACCACCGGACGAACGCATC </p>

*ermE**p

ATGGGGACCTCCTGGGGTGC GTTGGACCGCTGGATCCTACCAACCGGCACGATTGTGCCAC
AACAGCATCGCGGTGCCACGTGTGGACCGCGTCGGTCAGATCCTCCCCGCACCTCTCGCCAG
CCGTCAAGATCGAACCGCGTGACCCACATATGTTACCAATGCTTAATCAGTGAGGCACCTATCT
CAGCGATCTGTCTATTTTCGTTTCATCCATAGTTGCCTGACTCCCCGTCGTGTAGATAACTACGAT
ACGGGAGGGCTTACCATCTGGCCCCAGTGTGCAATGATACCGCGAGACCCACGCTCACCGG
CTCCAGATTTATCAGCAATAAACAGCCAGCCGGAAGGGCCGAGCGCAGAAGTGGTCCTGCA
ACTTTATCCGCCTCCATCCAGTCTATTAATTGTTGCCGGGAAGCTAGAGTAAGTAGTTCCGCCAG
TTAATAGTTTGCGCAACGTTGTTGCCATTGCTACAGGCATCGTGGTGTACGCTCGTCGTTTGG
TATGGCTTCATTACAGCTCCGTTCCCAACGATCAAGGCGAGTTACATGATCCCCCATGTTGTGC
AAAAAGCGGTTAGCTCCTTCGGTCCCTCCGATCGTTGTCAGAAGTAAGTTGGCCGCAGTGTTA
TCACTCATGGTTATGGCAGCACTGCATAATTCTCTTACTGTGCATGCCATCCGTAAGATGCTTTTC
TGTGACTGGTGAGTACTCAACCAAGTCATTCTGAGAATAGTGTATGCGGCGACCGAGTTGCTC
TTGCCCGGCGTCAATACGGGATAATACCGCGCCACATAGCAGAACCTTTAAAAAGTGTCTCATTT
GGAAAACGTTCTTCGGGGCGAAAACCTCTCAAGGATCTTACCCTGTTGAGATCCAGTTTCGATG
TAACCCACTCGTGCACCCAACTGATCTTCAGCATCTTTTACTTTACCAGCGTTTCTGGGTGAG
CAAAAACAGGAAGGCAAAATGCCGCAAAAAGGGAATAAGGGCGACACGGAAATGTTGAATAC
TCATACTCTTCCTTTTTCAATCATGATTGAAGCATTATCAGGGTTATTGTCTCATGAGCGGATA
CATATTTGAATGCTCGAG

Table S3, related to STAR Methods section “Bioactivity testing of streptophenazines and MIC determination.” Minimum Inhibitory Concentrations (MICs) ($\mu\text{g/mL}$) of oxo-streptophenazine A (9), streptophenazine C (13), streptophenazine A (16), and streptophenazine Q (20).

Strain	Compound 9	13	16	20
Group A <i>Streptococcus</i>	>50	>50	>50	2.5
<i>Acinetobacter baumannii</i> 5075	>50	>50	>50	40
<i>Klebsiella pneumoniae</i> 1100	>50	>50	>50	>40
MRSA TCH1516	>50	>50	>50	40

Table S4, related to STAR methods. Plasmids and strains used in this work.

Plasmids	Description	Sources
pCAP03	TAR cloning and broad-host-range heterologous expression vector; <i>CEN6-ARS4, oriT, traJ, pUC ori, Kan^r, Apra^r, pADH1, URA3, TRP1</i>	Tang et al., 2015
pSMM	Derivative of pCAP03 with 48 kb captured <i>spz</i> cluster	This work
pKDB01	Derivative of pCAP03 with 37.5 kb captured <i>spz</i> cluster	This work
pKDB02	Derivative of pKDB01 without TetR regulatory gene (<i>spz28</i>)	This work
pKDB01 Δ <i>spz3</i>	Derivative of pKDB01 without <i>spz3</i> (LysR type regulatory gene)	This work
pKDB02 Δ <i>spz3</i>	Derivative of pKDB02 without <i>spz3</i> (LysR type regulatory gene)	This work
pKDB01- <i>ermE</i> * <i>p-spz24</i>	Derivative of pKDB01 with <i>ermE</i> * promoter in front of <i>spz24</i> (LuxR-type regulatory gene), <i>Kan^r, Amp^r</i>	This work
pKDB03	Derivative of pKDB01 with <i>sp44-p21</i> cassette	This work
pKDB04	Derivative of pKDB03 with <i>act1p</i> cassette	This work
pKDB05	Derivative of pKDB04 with <i>ermE</i> * <i>p</i> cassette	This work
pKDB03 Δ <i>spz15</i>	Derivative of pKDB03 with <i>spz15</i> (putative adenylation protein-encoding gene) deleted	This work
pCAP03- <i>ermE</i> * <i>p</i>	Derivative of pCAP03 containing <i>ermE</i> * <i>p</i> cassette between XhoI and NdeI, <i>Kan^r, Amp^r</i>	This work
pCAP03- <i>act1p</i>	Derivative of pCAP03 containing <i>act1p</i> cassette between XhoI and NdeI, <i>Kan^r, Apra^r</i>	This work
pCAP03- <i>sp44-p21</i>	Derivative of pCAP03 containing <i>sp44-p21</i> cassette between XhoI, NdeI, <i>Kan^r, Apra^r</i>	This work
pIJ790	λ -RED (<i>gam, bet, exo</i>), <i>cat, araC</i> , rep101 ^{ts} , oriR101, P araBAD	Gust, 2003
pUB307	Self-transmissible plasmid that mobilizes other plasmids <i>in trans</i> for DNA transfer into hosts: RP4, <i>neo</i>	Flett, 1997
Strains	Description	
<i>Streptomyces</i>		
<i>Streptomyces</i> sp. CNB-091	Native producer of streptopenazines	Trischman et al., 1994
<i>Streptomyces coelicolor</i> M1146	Host strain for heterologous expression derived from <i>S. coelicolor</i> M145: Δ <i>act, Δred, Δcpk, Δcda</i> .	Gomez-Escribano and Bibb, 2011
<i>S. coelicolor</i> M1146-pCAP03	Heterologous host containing empty pCAP03 as a control	This work
<i>S. coelicolor</i> M1146-pSMM	Heterologous host containing pSMM (<i>spz</i> BGC captured in 48 kb DNA fragment)	This work
<i>S. coelicolor</i> M1146-pKDB01	Heterologous host with integrated pKDB01 (37.5 kb captured <i>spz</i> cluster)	This work
<i>S. coelicolor</i> M1146-pKDB02	Heterologous host with integrated pKDB02 (Δ <i>spz28</i>)	This work
<i>S. coelicolor</i> M1146-pKDB01 Δ <i>spz3</i>	Heterologous host with integrated pKDB01 Δ <i>spz3</i>	This work
<i>S. coelicolor</i> M1146-pKDB02 Δ <i>spz3</i>	Heterologous host with integrated pKDB02 Δ <i>spz3</i>	This work
<i>S. coelicolor</i> M1146-pKDB01- <i>ermE</i> * <i>p-spz24</i>	Heterologous host with integrated pKDB01- <i>ermE</i> * <i>p-spz24</i>	This work

<i>S. coelicolor</i> M1146-pKDB03	Heterologous host with integrated pKDB03	This work
<i>S. coelicolor</i> M1146-pKDB04	Heterologous host with integrated pKDB04	This work
<i>S. coelicolor</i> M1146-pKDB05	Heterologous host with integrated pKDB05	This work
<i>S. coelicolor</i> M1146-pKDB03 Δ spz15	Heterologous host with integrated pKDB03 Δ spz15	This work
<i>Escherichia coli</i>		
DH10B	F- <i>mcrA</i> Δ (<i>mrr-hsdRMS-mcrBC</i>), Φ 80 <i>lacZ</i> Δ M15, Δ <i>lacX74 recA1 endA1 araD139 Δ (<i>ara leu</i>)7697 <i>galU galK rpsL nupG</i> λ-. Storage and maintenance</i>	
BW25113	K-12 derivative: Δ <i>araBAD</i> , Δ <i>rhaBAD</i> ,	Datsenko and Wanner, 2000
BT340	DH5 α /pCP20, containing FLP recombinase	Cherepanov and Wackernagel, 1995
ET12567	F- <i>dam13</i> ::Tn9, <i>dcm6</i> , <i>hsdM</i> , <i>hsdR</i> , <i>recF</i> ,143 <i>zjj-202</i> ::Tn10, <i>galK2</i> , <i>galT22</i> , <i>ara-14</i> , <i>pacY1</i> , <i>xyl-5</i> , <i>leuB6</i> , <i>thi-1</i> , <i>tonA31</i> , <i>rpsL136</i> , <i>hisG4</i> , <i>tsx-78</i> , <i>mtl-1</i> , <i>glnV44</i> . Donor strain for conjugation between <i>E. coli</i> and <i>Streptomyces</i> in triparental mating	MacNeil et al., 1992
Other		
Group A <i>Streptococcus</i>	Clinical isolate used for bioactivity assays	Nizet lab, UCSD
<i>Acetivobacter baumannii</i> 5075	Clinical isolate used for bioactivity assays	Nizet lab, UCSD
<i>Klebsiella pneumoniae</i> 1100	Clinical isolate used for bioactivity assays	Nizet lab, UCSD
MRSA TCH1516	Clinical isolate used for bioactivity assays	Nizet lab, UCSD
<i>S. cerevisiae</i> VL6-48N	MAT α <i>trp1</i> - Δ 1 <i>ura3</i> - Δ 1 <i>ade2</i> -101 <i>his3</i> - Δ 200 <i>lys2 met14</i> <i>cir</i> ^o , TAR cloning	Kouprina and Larinova, 2016

Table S5, related to STAR methods. Primers used in this work.

Primer name	Sequence (5' to 3')	Purpose
Frag1_F	GAGTAGCAGCACGTTTCCTTATATGTAGCTTTCGACAT ATGCATGAGCTGTCTCCTGGTGTGGTGGGCAGG	TAR cloning <i>spz</i> BGC
Frag1_R	CGACCTGCCCCAACTCGACGGGCTGGAAGT	TAR cloning <i>spz</i> BGC
Frag2_F	CACGATGCCAGCAGCAGTCCCATGTCGTGG	TAR cloning <i>spz</i> BGC
Frag2_R	GCATCCAGTGCAGTACAGCCTCGCCGAGCG	TAR cloning <i>spz</i> BGC
Frag3_F	CCAGGTAATCCTCCAAGTCTCGCCGGTCCG	TAR cloning <i>spz</i> BGC
Frag3_R	CTCATAAGGATGCCTTCTGCGGGTGTGGAGACC	TAR cloning <i>spz</i> BGC
Frag4_F	CGATCGCCGCCTGGACCGACTGGAGCAGG	TAR cloning <i>spz</i> BGC
Frag4_R	GCTGCCGTGAACGCCAACGGGAAGGTGGAC	TAR cloning <i>spz</i> BGC
Frag5_F	GAGGGCGGTCAGGGCAGGATTCACGGAATGC	TAR cloning <i>spz</i> BGC
Frag5_R	GGAGGCAGGCGCTACTGGACGTTCTCATCC	TAR cloning <i>spz</i> BGC
Frag6_F	CCATCATCTTGAGCTCTTCGTGACGACCATG	TAR cloning <i>spz</i> BGC
Frag6_R	GAACCAGGAGATCGCCGACCGGCTCGTCCT	TAR cloning <i>spz</i> BGC
Frag7_F	CTCGGCCTGGAATCACTTCAGATGATGCGCCTG	TAR cloning <i>spz</i> BGC
Frag7_R	GTTCTGATCGCCAGTTCCTGCAGAGCGTG	TAR cloning <i>spz</i> BGC
Frag8_F	GTCACTCCGACGAACATCAGGATGGAGGCCG	TAR cloning <i>spz</i> BGC
Frag8_R+TetR	CTCGGTTTGACGCCTCCCATGGTATAAATAGTGGCTC GAGGCACAGGACGGCATCGCCCGAGCTGAGC	TAR cloning <i>spz</i> BGC
Frag8_R-TetR	CTCGGTTTGACGCCTCCCATGGTATAAATAGTGGCTC GAGGGTCCGTCCGGTCCCAGGAGAACCAGCTGATC GGCTTGCGGAGAGACCGGGCGGGGATCAGCCGAGG AGGGCGGGCGTGTACGCCAATCGACTGGCGAGC GGCAT	TAR cloning <i>spz</i> BGC
Δ<i>spz3</i>_F	GCTCGTCGCACCGTCGCCCGCTCCTGAACCGGCC GCGGGGAAGGACCACGGTTCATGTGCAGCTCCATCA GCAAAAG	Delete <i>spz3</i> (LysR)
Δ<i>spz3</i>_R	GTCGTCGCACCGTCGCCCGCTCCTGAACCGGCC GCGGGGAAGGACCACGGTTCATGTGCAGCTCCATCA GCAAAAG	Delete <i>spz3</i> (LysR)
Δ<i>spz3</i>_seq_F	GTGAAGGCGACGAACAGGAAGTG	Confirm <i>spz3</i> deletion
Δ<i>spz3</i>_seq_R	GATCGAGCCC GCCGACGTG	Confirm <i>spz3</i> deletion
Upreg-<i>spz24</i>_F	GTCGGGAGCGTGCGGCACAGTGCGCACGTGGGCA GGACTGGCACGGTCATATGGGACCTCCTGGGGTG CGTTGG	Insert <i>ermE</i> * promoter in front of <i>spz24</i> (LuxR)
Upreg-<i>spz24</i>_R	CCCGATCGTGTACGTGCGCGCCGTGCGCCCTAGG CCCTCCGGGCGGACGCTCGAGCATTCAAATATGTAT CCGCTCATG	Insert <i>ermE</i> * promoter in front of <i>spz24</i> (LuxR)
Upreg-<i>spz24</i>_seq_F	CAGCAGGCTGCTCTTGCCGACAC	Sequencing to confirm insertion
Upreg-<i>spz24</i>_seq_R	GATGCGCGGTATGTGGGAGCGC	Sequencing to confirm insertion
<i>aac(3)IV</i> +FRT_F	ATTATACATATGTGCAGTTCGAAGTTCCTATTCTCTAG AAAGTATAGGAATTCGGTTCATGTGCAGCTCCATCA GCAAAAG	Amplify <i>aac(3)IV</i> gene with FRT sites for cloning between pET28a NdeI and HindIII sites
<i>aac(3)IV</i> +FRT_R	ATTATAAAGCTTGAAGTTCCTATACTTTCTAGAGAATA GAACTTCTCAGCCAATCGACTGGCGAGCGGCATCG	Amplify <i>aac(3)IV</i> gene with FRT sites for cloning between pET28a NdeI and HindIII sites
<i>sp44</i>_F	GACGCCTCCCATGGTATAAATAGTGGCTCGAGGGTG AACCGATCTCCTCGTTAGGGTC	Amplify <i>sp44</i> promoter with homology to pCAP03 and <i>aac(3)IV</i>
<i>sp44</i>_R	TCTAGAGAATAGGAACTTCGAAGTGCACATATGATAG TGGCTGTTACATTGCAACCGTCTCTG	Amplify <i>sp44</i> promoter with homology to pCAP03 and <i>aac(3)IV</i>

<i>p21_F</i>	CTATTCTCTAGAAAGTATAGGAACTTCAAGCTTTGCT CGAGTGTGCGGGCTCTAACACGTC	Amplify <i>p21</i> promoter with homology to pCAP03 and <i>aac(3)IV</i>
<i>p21_R</i>	AGCACGTTCCCTTATATGTAGCTTTTCGACATATGGTCC GTACCTCCGTTGCTCGACTAGAT	Amplify <i>p21</i> promoter with homology to pCAP03 and <i>aac(3)IV</i>
<i>aac(3)IV+FRT_sp44_pCAP03_F</i>	CAGAGACGGTTTGAATGTGAACAGCCACTATCATAT GTGCAGTTCGAAGTTCTATTCTCTAGA	Amplify <i>aac(3)IV</i> +FRT with homology to <i>sp44</i> promoters
<i>aac(3)IV+FRT_p21_pCAP03_R</i>	ACGTGTTAGAGCCCCGACACTCGAGCAAAGCTTGAA GTTCTTACTACTTTCTAGAGAATA	Amplify <i>aac(3)IV</i> +FRT with homology to <i>p21</i> promoter
<i>sp44+p21_cas_insert_F</i>	TGGTCGGCGTCTCTGTTGCTGAAACTGATGCTCTGG GCGACGTGCGGCATCTCGAGGGTGAACCGATCTCCT CGTTAGG	Amplify <i>sp44-p21</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>sp44+p21_cas_insert_R</i>	CCGGGGACGTGGGTGCCGAGTCCGGCGAGGTAGGT GTTCTCGAACTTCATCATATGGTCCGTACCTCCGTTG CTCGAC	Amplify <i>sp44-p21</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>sp44+p21_cas_insert_seq_F</i>	GTGAACAGGAGATGCCGGGTG	Confirm <i>sp44-p21</i> cassette insertion
<i>sp44+p21_cas_insert_seq_R</i>	GCCTCGTACCAGCCTTCCTG	Confirm <i>sp44-p21</i> cassette insertion
<i>actlp_F</i>	CCGGTCCGTGAACGCGGTGGAGCCCTATGTCTCTTA AGTGTCCCTCCGTGCC	Amplify <i>actI</i> promoter with homology to <i>actII</i> -ORF4 and pCAP03
<i>actlp_R</i>	CACGTTCCCTTATATGACTTTTCGAGATGCGTTCGTC CGGTGGTCGTCGTC	Amplify <i>actI</i> promoter with homology to <i>actII</i> -ORF4 and pCAP03
<i>actII-ORF4_F</i>	CCGCTCGCCAGTCGATTGGCTGACATATGCCACTGC CTCTCGGTAATAATCC	Amplify <i>actII</i> -ORF4 with homology arms for assembly with <i>aac(3)IV</i> and <i>act1p</i>
<i>actII-ORF4_R</i>	GGCAGGGAGGGGAACACTTAAGAGACATAGGGCTC CACCGCGTTCACGGACCGG	Amplify <i>actII</i> -ORF4 with homology arms for assembly with <i>aac(3)IV</i> and <i>act1p</i>
<i>aac(3)IV+actlp_F</i>	GACGCCTCCCATGGTATAAATAGTGGCTCGAGGGTT CATGTGCAGCTCCATCAGCAAAAG	Amplify <i>aac(3)IV</i> with homology to <i>actII</i> -ORF4 and pCAP03
<i>aac(3)IV+actlp_R</i>	GGATTTTACCGAGAGGCAGTGGCATATGTCAGCCAA TCGACTGGCGAGCGG	Amplify <i>aac(3)IV</i> with homology to <i>actII</i> -ORF4 and pCAP03
<i>actlp_cas_insert_F</i>	GGGTCCGCGTGTACGGAAGGTCAAGAATCTTCGGG TCGGCGGAAGCCACGATGCGTTCGTCGGTGGTTCG TCGTC	Amplify <i>actlp</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>actlp_cas_insert_R</i>	GTGCTCGAATGTCCATACACCCAAGACGTAGAAGT TCTCTGGAGGAACGACTCGAGGGTTCATGTGCAGCT CCATCAG	Amplify of <i>actlp</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>actlp_cas_insert_seq_F</i>	CTCGAGCCAGTAGCGGGATC	Confirm <i>actlp</i> cassette insertion within <i>spz</i> BGC
<i>actlp_cas_insert_seq_R</i>	CATTGTGCACGGTCCACCG	Confirm <i>actlp</i> cassette insertion within <i>spz</i> BGC
<i>ermE*p_F</i>	GTGCCTCACTGATTAAGCATTGGTAACATATGTGGTG CACGCGGTGATCTTGACGGCTG	Amplify <i>ermE*</i> promoter with homology arms for assembly with <i>bla</i> and pCAP03
<i>ermE*p_R</i>	CAGCACGTTCCCTTATATGTAGCTTTTGAATGGGGACC TCCTGGGGTGCGTTGGACC	Amplify <i>ermE*</i> promoter with homology arms for assembly with <i>bla</i> and pCAP03
<i>bla+ermE*p_F</i>	TGACGCCTCCCATGGTATAAATAGTGGCTCGAGCATT CAAATATGTATCCGCTCATGAGA	Amplify <i>bla</i> with homology to <i>ermE*</i> and pCAP03 for assembly
<i>bla+ermE*p_R</i>	CAGCCGTCAAGATCGACCGCGTGCACCACATATGTT ACCAATGCTTAATCAGTGAGGCAC	Amplify <i>bla</i> with homology to <i>ermE*</i> and pCAP03 for assembly
<i>ermE*p_cas_insert_F</i>	CAGGGTCGCCGTGACGGCGAACTGGAGGCGGGCGA GGTCGACGACGTCCACATGGGGACCTCCTGGGGTG CGTTGG	Amplify <i>ermE*p</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>ermE*p_cas_insert_R</i>	CTTCCACAGCCACGCCCGCGCCCTCCTCGGCTGATC CCCGCCCGGTCTCTCTCGAGCATTCAAATATGTATCC GCTCATG	Amplify <i>ermE*p</i> cassette with homology sequences for targeted insertion into <i>spz</i> BGC
<i>ermE*p_cas_insert_seq_F</i>	GTGAGGTGCGCTCGTGCGCAG	Confirm <i>ermE*p</i> cassette insertion within <i>spz</i> BGC
<i>ermE*p_cas_insert_seq_R</i>	GTCCTGCCCGAGCAGTATGTCC	Confirm <i>ermE*p</i> cassette insertion within <i>spz</i> BGC
Δ <i>spz15_F</i>	ACGGGTTCCAGGCGGCCGGAGACGACGTCCTCGCG CAGCGCCCCCGGTCTCAGCCAATCGACTGGCGAG CGGCAT	Delete <i>spz15</i> (adenylation protein)
Δ <i>spz15_R</i>	GGCGGGTCCGTCGGGCTCCAGCCCTGAGCCCCGGGGC ACTCCTCCGCCTCCGGTTCATGTGCAGCTCCATCAG CAAAAG	Delete <i>spz15</i> (adenylation protein)
Δ <i>spz15_seq_F</i>	GTCAGGGCAGGATTCACGGAATG	Confirmation of <i>spz15</i> deletion
Δ <i>spz15_seq_R</i>	GTACGAACTCCGTTTCGACAGGGTG	Confirmation of <i>spz15</i> deletion

Table S6, related to STAR methods. Primers used for RT-PCR.

Primer name	Sequence (5' to 3')	Primer name	Sequence (5' to 3')
<i>hrdB_F</i>	CCTCCGCCTGGTGGTCTCG	<i>hrdB_R</i>	AACTTGTAGCCCTTGGTGTAGTCGAAC
<i>spz1_F</i>	CGTCACCGAGAGCAGCCACAGCG	<i>spz1_R</i>	GTTGCTCCAAGCACTGCGACTGCC
<i>spz2_F</i>	CGATCAGGGAGAAGTTCACGACG	<i>spz2_R</i>	CTGCACTTCTGTTCGTGCGCCTT
<i>spz3_F</i>	GGACATACTGCTCGGGCAGGAC	<i>spz3_R</i>	CAGTTGGAGTACTTCTCGCGG
<i>spz4_F</i>	GCCTGATCAGCACCACGTCGG	<i>spz4_R</i>	GAAGTTCGCGTACGACCTGTGCG
<i>spz5_F</i>	CGATCTGCTCGAAGTCGAACAGCC	<i>spz5_R</i>	GACACCATGTACGCGTCACTCATCC
<i>spz6_F</i>	GTGCAGAGTCTCCAGAACGTCGTGG	<i>spz6_R</i>	CACACCGCAGGTCGTGTTACCTTC
<i>spz7_F</i>	GCGAGGCTGTACTCGCACTGGATG	<i>spz7_R</i>	GTCTCCGAATACGCGCTGGGCAC
<i>spz8_F</i>	GAACGTGGTCACCGAATGTCCAG	<i>spz8_R</i>	GACCGAGCTCGACCTGTACTGCTTC
<i>spz9_F</i>	CTCGAAGAAGCTCTCGTCCGGTTCG	<i>spz9_R</i>	CTTGAGGACTGGAAGCCGTACGGG
<i>spz10_F</i>	GTGGTCTGAACAGGTCGATCAGCTC	<i>spz10_R</i>	GAGAAGGATTCCTGACCCTCGTCCG
<i>spz11_F</i>	GTCTTGCTCCGCATCAGGAATTCC	<i>spz11_R</i>	GATCACCTCGACATGATCGTCGAC
<i>spz12_F</i>	CTCACGTGGTGCATGACCACCTGG	<i>spz12_R</i>	GACATCGTCGGCGTCTTCGAGAG
<i>spz13_F</i>	CGATGATCGCCTGGTGGGATCTC	<i>spz13_R</i>	GTGGTGGCGACGACGTCTGAAC
<i>spz14_F</i>	CGAGATGACGTACGTCCACCACGC	<i>spz14_R</i>	GCTTCCAGACGGCCTACCGGCATC
<i>spz15_F</i>	CACGTCCGCCAGGTGGTAGGTGAG	<i>spz15_R</i>	GAGGTGACGGTGTGGTTCTCGGTG
<i>spz16_F</i>	CTCCACCCTGTCTGAACCGGAGTTC	<i>spz16_R</i>	GACCGCCGATATCACCGTCCCCTTTC
<i>spz17_F</i>	GCGAGAACATCCTGCTGCGCCAC	<i>spz17_R</i>	CTCCTATGTGGTGTGTCGACGCCTTC
<i>spz18_F</i>	CGGTTGATCTCGACGACACCGAGC	<i>spz18_R</i>	CTGGACTCCGCGATCCTCATCCGC
<i>spz19_F</i>	CACCAGGAAGGTGCGGATGTCTGTG	<i>spz19_R</i>	CTGGAATCCCCGATCCACCCG
<i>spz20_F</i>	CAGGTACTTGTCTCCGACGCCCGAC	<i>spz20_R</i>	CTCGACTACGAGCTGCCGATGCTG
<i>spz21_F</i>	CTGGTCTGAACGGGTTTCATGAAGTCG	<i>spz21_R</i>	GCTGACCCGGCATCTCTGTTTAC
<i>spz122_F</i>	CTGGTCAACCGGAGAGGCGGTC	<i>spz22_R</i>	GATTCCAGATCGAAGCCGATCGGATC
<i>spz23_F</i>	CTAGTCGCCGGAGAGGAGCTGCGC	<i>spz23_R</i>	GTATCGCGTCTCTGCTCCTCAGC
<i>spz24_F</i>	CAACTCCTCGTAGGCGTGCGACAC	<i>spz24_R</i>	CTACCGTCTCTCCTGCGGCTGGATC
<i>spz25_F</i>	GGTGGTGTAGGTGAAGTCCACCCAC	<i>spz25_R</i>	GAAGTCGTACGACCAGAGCGGTGTG
<i>spz26_F</i>	CGAACATCAGGATGGAGGCCGTCAG	<i>spz26_R</i>	GCTGGAGTACTTCTGGTGGGGCTCG
<i>spz27_F</i>	CGTCGAGCCGGAGGTCTGAAGAAG	<i>spz27_R</i>	CAAGGACATCTGGCCGAGTGGCTG
<i>spz28_F</i>	CGATCAGCTGGTCTCTCTGGGACCG	<i>spz28_R</i>	GACTAGGGCCAGGCCAAGGAGCG